

MCI-Gel[®] ProtEx, CQA and CQK ion-exchange columns

Chromatograms illustrating the performance of the MCI-Gel ProtEx ion-exchange columns, followed by some chromatograms of the MCI-Gel CQA and CQK columns. Type of column and sample are shown in each graph.

Conditions

Column : MCI GEL® ProtEx-DEAE 4.6mm I.D.×50mm

Eluent : A 20mM Tris-HCl pH8.0

B A+0.5M NaCl

A → B 30min linear gradient

Flow rate : 0.5 ml/min


Column temp. : ambient

Detection : 280nm

Sample : 1. Myoglobin 25μg

2. Conalbumin 25μg

3. Trypsin Inhibitor 50μg



Conditions

Column : MCI GEL® ProtEx-DEAE 4.6mm I.D.×50mm

Eluent : A 20mM Tris-HCl pH8.0

B A+0.5M NaCl

A → 10% B 30min linear gradient

Flow rate : 0.5 ml/min


Column temp. : ambient

Detection : 280nm

Sample : 1. Hb A₂ 100μg

2. Hb S 100μg

3. Hb A₀ 100μg



Conditions

Column : MCI GEL® ProtEx-DEAE 4.6mm I.D.×50mm

Eluent : A 20mM Tris-HCl pH8.0

B A+0.5M NaCl


: 5% B → 70% B 30min linear gradient

Flow rate : 0.5 ml/min

Column temp.: ambient

Detection : 280nm

Sample : recombinant hGH 10µg



Conditions

Column : MCI GEL® ProtEx-DEAE 4.6mm I.D.×50mm

Eluent : A 20mM Tris-HCl pH8.15

B A+0.5M NaCl


A → 50% B 30min linear gradient

Flow rate : 0.5 ml/min


Column temp. : ambient

Detection : 280nm

Sample : recombinant epidermal growth factor (EGF)





Conditions
 Column : MCI GEL® ProtEx-DEAE 4.6mm I.D.×50mm
 Eluent : A 20mM Trimethylenediamine-HCl pH9.75
 B A+0.5M NaCl
 Flow rate : 0.5 ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : recombinant IL-2 1.5µg
 1. BSA (stabilizer) 400µg
 2. IL-2 (Met-ox)
 3. IL-2



Conditions

Column : MCI GEL® ProtEx-SP 4.6mm I.D.×50mm
 Eluent : A 20mM Phosphate buffer pH6.0
 B A+0.5M NaCl
 A → B 20min linear gradient
 Flow rate : 0.5 ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : 1. Ribonuclease A 10µg
 2. α-Chymotrypsinogen A 5µg
 3. Cytochrome C 5µg





Conditions

Column : MCI GEL® ProtEx-SP 4.6mm I.D.×50mm

Eluent : A 20mM Bis-Tris HCl buffer pH6.0

B A+0.5M NaCl

7% B → 40% B 20min linear gradient

Flow rate : 0.5 ml/min


Column temp. : ambient

Detection : 415nm

Sample : GHb

1. GHb A1c

2. Hb A0



Conditions

Column : MCI GEL® ProtEx-DEAE 4.6mm I.D.×50mm

Eluent : A 20mM HEPES buffer pH7.6

B A+0.5M NaCl

A100% → B45% in 30min. B45% for 5min


B45% → B100% in 5min. B100% for 10min

Flow rate : 0.5 ml/min

Column temp. : 25°C

Detection : 280nm

Sample : IgG2b, κ (mouse) 10 μ g



Conditions

Column : MCI GEL® ProtEx-DEAE 4.6mm I.D.×50mm

Eluent : A 10mM HEPES buffer pH8.0

B A+0.5M NaCl


A100% → B100% in 30min. linear gradient

Flow rate : 0.5 ml/min

Column temp. : 25°C

Detection : 280nm

Sample : IgG1 MOPC21 (mouse) 10µg



Conditions
 Column : MCI GEL® CQA31S 7.5mm I.D.×75mm
 MCI GEL® CQA35S 7.5mm I.D.×75mm
 Eluent : A 14mM Tris-HCl buffer pH8.2
 B A +0.5M NaCl
 A → B 30min linear gradient
 Flow rate : 1.0 ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : 1. Myoglobin 60µg
 2. Ovalbumin 200µg
 3. Trypsin Inhibitor 200µg

Conditions
 Column : MCI GEL® CQA31S 7.5mm I.D.×75mm
 MCI GEL® CQA35S 7.5mm I.D.×75mm
 Eluent : A 14mM Tris-HCl buffer pH8.2
 B A +0.5M NaCl
 A → B 30min linear gradient
 Flow rate : 1.0 ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : 1. Myoglobin 120µg
 2. Transferrin 160µg
 3. β-Lactoglobulin 400µg








Fig. 4-20 Separation of protein mixture



Conditions
 Column : MCI GEL® CQK30S 7.5mm I.D.×75mm
 MCI GEL® CQK31S 7.5mm I.D.×75mm
 Eluent : A 20mM Phosphate buffer pH6.8
 B A +0.5M NaCl
 A → B 30min linear gradient
 Flow rate : 1.0 ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : 1. Myoglobin 50µg
 2. Ribonuclease A 200µg
 3. Cytochrome C 60µg
 4. Lysozyme 80µg




Conditions
 Column : MCI GEL® CQK30S 7.5mm I.D.×75mm
 MCI GEL® CQK31S 7.5mm I.D.×75mm
 Eluent : A 20mM Phosphate buffer pH6.8
 B A +0.5M NaCl
 A → B 30min linear gradient
 Flow rate : 1.0 ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : 1. Trypsinogen 80µg
 2. Ribonuclease A 90µg
 3. α-Chymotrypsinogen A 60µg

Fig. 4-22 Separation of lipoxidase

Conditions
Column : MCI GEL® CQA31S 7.5mm I.D.×75mm
Eluent : A 14mM Tris-HCl buffer pH8.2
 B A +0.5M NaCl
 A → B 30min linear gradient
Flow rate : 1.0 ml/min
Column temp. : ambient
Detection : 280nm
Sample : Lipoxidase

Product from company A



Product from company B

